silasoft® Normal

Instructions for use

**Precision impression material, C-silicone based,
low consistency – light bodied**

Fig. 1

1. Preparation of the impression

Depending on the impression technique select an appropriate tray. For optimal adhesion of the impression to the tray apply a thin layer of **Sili Adhesive Spray** or **Sili Adhesive liquid** onto the tray and let dry for 60 sec (Fig. 1). Take preliminary impression with **silaplast FUTUR**. To facilitate the flowing out of excessive light body material from the impression cutting of escape vents and removal of interdental septi in the impression with the use of the special instrument **Deta-Cut** is recommended.


Fig. 2

2. Mixing

Extrude **silasoft® Normal** and **cat[®] universal FUTUR** in equal lengths and thickness onto mixing pad (Fig. 2). Using catalyst liquid, the mixing ratio is 1 cm of **silasoft® N** to 1 drop of catalyst. Mix according to dosing ratio for approx. 30 sec. with spatula until homogeneous uniform colouring is achieved. Load preliminary impression. A normal syringe or disposal syringe can be used for bubble free application (Fig. 3 + 4).

Fig. 3

3. Disinfection

After removal rinse impression under lukewarm water. Disinfection for 15 min. in 2 % glutaraldehyde.


Fig. 4

4. Model casting

Casting between 30 minutes and 24 hours after removal of the impression. Recommended materials: Dental plasters of class III and IV-also commercially available model acrylics.

5. Galvanisation

Impressions can be electroplated with commercially available copper and silver plating baths.


6. Cleaning of trays

Remove set material with a blunt instrument. **Sili Spray** or **Sili liquid** film can be easily removed with **Sili Spray Solvent** (free of FCHC). Work in well ventilated rooms. Clean und disinfect as usual.

Important working hints

- **silasoft® N Bottle: Mix contents homogeneously prior to the first use.**
- Increased temperatures accelerate, decreased temperatures retard the setting time.
- Catalyst bottles or tubes should be closed tightly immediately after use.
- Don't leave any residual material in the patient's mouth.
- Avoid staining of clothing. Hardened impression materials cannot be removed.
- Avoid contact of catalyst with eyes, may cause irritation. If necessary flush eyes with plenty of water and seek immediate medical care.

silasoft® N / cat[®] universal FUTUR:

Indications for use:

- Corrective ("Wash")-material for double impression
- Reline impression

Technical Data:

DIN EN ISO 4823 – Type 3

- **Mixing ratio:**
silasoft® N / cat[®] universal FUTUR
1 cm silasoft® N = 0.50 g /
1 cm cat[®] = 0.05 g
silasoft® N / silasoft® cat[®]
1 cm silasoft® N = 0.50 g /
1 drop cat[®] = 0.02 g
- **Colour code:**
silasoft® N: rose-opaque
cat[®] universal FUTUR: blue
silasoft® cat[®]: light blue
- **Mixing time:** approx. 30 sec.
- **Working time:** 1 min. 15 sec*
- **Setting time:** approx. 3 min. 45 sec.*
- **Time in mouth:** approx. 2 min. 30 sec.
- **Strain in compression:** 11.20 %
- **Recovery from deformation:** 98.40 %
- **Linear dimensional change:** 0.50 %
- **Application:** At 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5 % rel. humidity

* from beginning of mixing at 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5 % rel. humidity. Increased temperatures accelerate, decreased temperatures retard a. m. times.

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