




The Inventors of Hardened Plastic Teeth.



JUSTI BLEND®

MOLD CHART

ORTHOCHROMATIC ANTERIOR &
INTERCUSPAL-OCCLUSAL DESIGNED POSTERIOR
MULTI-LAYERED
HARDENED PLASTIC TEETH























for Today's Cosmetic Industry



Highly Cross linked/Intersectional Polymer Network • Wear Resistant • Vacuum Processed • Fluorescent/Photogenic
No diatoric retention necessary. 100% bondability with methyl methacrylate denture base material.
Justi Blend complies with ANSI/ADA Specification #15, ISO 3336, BS 3990, NF S91-213, DIN 13914

JUSTI BLEND®

MULTI-BLENDED PLASTIC ANTERIORS
27 UPPERS & 8 LOWERS

| | |
|---|--|
| A4 10.0 8.3  49.0 43.78 4 | D6 12.0 9.0  54.0 47.68 6 |
| A7 10.0 8.8  52.5 46.98 7 | E1 9.3 7.3  45.0 39.86 1 |
| B3 9.5 8.0  49.0 42.30 3 | E2 10.0 7.8  46.5 41.31 2 |
| B8 11.3 8.8  54.0 47.09 8 | E3 10.5 8.3  50.0 43.29 3 |
| C1 10.0 7.5  45.0 40.30 1 | E6 10.0 8.5  52.0 46.01 6 |
| C2 10.0 8.0  47.5 41.75 2 | E8 11.0 8.8  55.0 48.90 8 |
| C5 11.8 8.8  51.0 45.51 5 | F1 8.5 7.8  46.0 39.29 1 |
| C6 10.8 8.5  51.0 43.76 6 | F4 9.0 8.5  50.0 43.75 4 |
| D1 10.0 7.5  45.0 38.36 1 | G5 10.5 8.8  51.5 44.78 5 |
| D2 10.0 8.0  47.0 42.24 2 | G7 10.8 8.8  54.0 47.43 7 |
| D3 10.5 8.3  49.0 43.44 3 | H2 10.0 8.0  47.5 41.06 2 |

Photographs of teeth are within plus or minus 3% of actual size.

**MULTI-BLENDED PLASTIC ANTERIORS
27 UPPERS & 8 LOWER**



| | | | |
|------------|--|--|---------------------------------|
| H4 | | | 4 |
| 10.5 8.3 | | | 51.0 44.67 () |
| J2 | | | 2 |
| 10.0 7.8 | | | 46.5 40.75 () |
| J3 | | | 3 |
| 10.0 8.0 | | | 48.0 42.68 () |
| J5 | | | 5 |
| 11.5 9.0 | | | 53.0 47.14 () |
| K4 | | | 4 |
| 10.8 8.0 | | | 49.5 43.90 () |
| I | | | CI, DI EI, FI |
| 9.0 4.8 | | | 34.5 30.29 () |
| 2 | | | C2, D2 E2, H2 J2 |
| 8.8 4.5 | | | 36.0 29.86 () |
| 3 | | | B3, D3 E3, J3 |
| 8.3 4.8 | | | 37.0 30.17 () |
| 4 | | | A4, F4 H4, K4 |
| 9.0 4.8 | | | 38.0 31.07 () |
| 5 | | | C5, G5 J5 |
| 9.8 5.3 | | | 40.0 32.80 () |
| 6 | | | C6, D6 E6 |
| 8.5 5.3 | | | 39.5 33.32 () |

| | | | |
|-----------|--|--|-------------------|
| 7 | | | A7, G7 |
| 9.5 5.0 | | | 41.0 33.38 () |

| | | | |
|-----------|--|--|-------------------|
| 8 | | | B8, E8 |
| 9.3 5.5 | | | 41.5 35.77 () |

0° POSTERIORS - 3 UPPERS & LOWER

| | | | |
|--------------|--|--|------|
| 29 MF | | | 6.69 |
| 28.5 | | | 8.15 |
| 30.5 | | | |

| | | | |
|--------------|--|--|------|
| 31 MF | | | 6.60 |
| 30.0 | | | 7.62 |
| 32.8 | | | |

| | | | |
|--------------|--|--|------|
| 33 MF | | | 7.31 |
| 32.2 | | | 8.30 |
| 34.0 | | | |

10° POSTERIORS - 3 UPPERS & LOWER

| | | | |
|-------------|--|--|------|
| 30 F | | | 8.41 |
| 30.0 | | | 9.04 |
| 33.0 | | | |

| | | | |
|-------------|--|--|------|
| 32 F | | | 8.72 |
| 31.0 | | | 9.63 |
| 33.5 | | | |

| | | | |
|-------------|--|--|-------|
| 34 F | | | 9.15 |
| 33.0 | | | 10.17 |
| 36.5 | | | |

**20° POSTERIOURS
6 UPPERS & LOWERES**

**COMPLIES WITH
ADA SPEC. #15**

**ANATOMICAL MEDIUM - 33°
3 UPPERS & LOWERES**

| | | | |
|-----------------|------|--|------|
| 29 L | | | 7.94 |
| | | | 8.51 |
| | 29.1 | | 31.5 |

| | | | |
|-----------------|------|--|------|
| 30 M | | | 8.92 |
| | | | 9.33 |
| | 29.9 | | 32.9 |

| | | | |
|-----------------|------|--|------|
| 31 L | | | 8.64 |
| | | | 8.75 |
| | 31.1 | | 33.5 |

| | | | |
|-----------------|------|--|------|
| 32 M | | | 8.88 |
| | | | 9.47 |
| | 31.9 | | 35.9 |

| | | | |
|-----------------|------|--|------|
| 33 L | | | 9.39 |
| | | | 9.72 |
| | 33.1 | | 35.5 |

| | | | |
|-----------------|------|--|-------|
| 34 M | | | 10.02 |
| | | | 10.67 |
| | 34.0 | | 38.0 |

| | | | |
|-----------------|------|--|------|
| 29 M | | | 7.66 |
| | | | 8.49 |
| | 29.1 | | 31.5 |

**ANATOMICAL LONG - 33°
3 UPPERS & LOWERES**



| | | | |
|-----------------|------|--|------|
| 30 L | | | 8.82 |
| | | | 9.20 |
| | 29.2 | | 31.9 |

| | | | |
|-----------------|------|--|------|
| 31 M | | | 8.10 |
| | | | 9.36 |
| | 31.1 | | 33.5 |

| | | | |
|-----------------|------|--|------|
| 32 L | | | 9.06 |
| | | | 9.58 |
| | 31.2 | | 34.2 |

| | | | |
|-----------------|------|--|------|
| 33 M | | | 8.38 |
| | | | 9.58 |
| | 33.1 | | 35.5 |

| | | | |
|-----------------|------|--|-------|
| 34 L | | | 9.68 |
| | | | 10.33 |
| | 33.6 | | 36.9 |



Mold Classification Symbols:
 ■ **SQUARE** ▼ **TAPERING** ● **OVOID**
 Combined Symbol: ▼● **Tapering Ovoid Type**



ORDER FORM

Order No.: _____ Date: _____

Bill To: _____

Address: _____

City: _____ State: _____ Zip: _____

Phone # _____

Ship To: _____

Address: _____

City: _____ State: _____ Zip: _____

1 x 6 27 UPPER ANTERIOR MOLDS Boxes

| Molds | 200 | 202 | 204 | 206 | 208 | 210 | 212 | 214 | 216 |
|-------|-----|-----|-----|-----|-----|-----|-----|-----|-----|
| A4 | | | | | | | | | |
| A7 | | | | | | | | | |
| B3 | | | | | | | | | |
| B8 | | | | | | | * | * | * |
| C1 | | | | | | | * | * | * |
| C2 | | | | | | | | | |
| C5 | | | | | | | * | * | * |
| C6 | | | | | | | | | |
| D1 | | | | | | | | | |
| D2 | | | | | | | | | |
| D3 | | | | | | | | | |
| D6 | | | | | | | * | * | * |
| E1 | | | | | | | * | * | * |
| E2 | | | | | | | * | * | * |
| E3 | | | | | | | | | |
| E6 | | | | | | | | | |
| E8 | | | | | | | * | * | * |
| F1 | | | | | | | * | * | * |
| F4 | | | | | | | | | |
| G5 | | | | | | | * | * | * |
| G7 | | | | | | | | | |
| H2 | | | | | | | | | |
| H4 | | | | | | | * | * | * |
| J2 | | | | | | | * | * | * |
| J3 | | | | | | | | | |
| J5 | | | | | | | * | * | * |
| K4 | | | | | | | | | |
| Total | | | | | | | | | |

1 x 6 8 LOWER ANTERIOR MOLDS Boxes

| Molds | 200 | 202 | 204 | 206 | 208 | 210 | 212 | 214 | 216 |
|-------|-----|-----|-----|-----|-----|-----|-----|-----|-----|
| 1 | | | | | | | | | |
| 2 | | | | | | | | | |
| 3 | | | | | | | | | |
| 4 | | | | | | | | | |
| 5 | | | | | | | * | * | * |
| 6 | | | | | | | | | |
| 7 | | | | | | | | | |
| 8 | | | | | | | * | * | * |
| Total | | | | | | | | | |

* Not Available
Anterior Cards (1 x 6) - 16 cards per box.

Posterior Cards (1 x 8) - 20 cards per box.

| | |
|-----------------------------|--|
| Total Upper Anteriors 1 x 6 | |
| Total Lower Anteriors 1 x 6 | |
| Total Anteriors 1 x 6 | |

| | |
|------------------------------|--|
| Total Upper Posteriors 1 x 8 | |
| Total Lower Posteriors 1 x 8 | |
| Total Posteriors 1 x 8 | |

Item #079-58-0001

Item #079-58-0003

Remarks _____

Signed _____

JUSTI BLEND PRODUCTS:

| Item No. | Description |
|---------------|---|
| 079-58-0001 | Justi Blend Anterior Teeth (1x6) |
| 079-58-0003 | Justi Blend Posterior Teeth (1x8) |
| 084-58-0004 | Justi Blend Combination (1x28) |
| 084-59-0026/C | Justi Blend Maxillary Centrals Shade Guide |
| 084-59-0026 | Justi Blend Extended Blend Selector in 9 shades with one Kroma Tab |
| 084-59-0010 | Justi Blend Anterior Working Mold Guide Contains one set of each mold - 35 sets, all one shade |
| 084-59-0011 | Justi Blend Posterior Working Mold Guide Contains one set of each mold - 36 posteriors, all one shade |
| 086-55-0101 | Tooth Cabinet - 10 drawer capacity - 210 sets of Justi Blend (anterior/posterior) |



Mold charts, rulers and prescription pads are also available.

1 x 8 18 UPPER POSTERIOR MOLDS Boxes

| Molds | 200 | 202 | 204 | 206 | 208 | 210 | 212 | 214 | 216 |
|-------------------------|-----|-----|-----|-----|-----|-----|-----|-----|-----|
| 0° | | | | | | | | | |
| 29MF | | | | | | | | | |
| 31MF | | | | | | | | | |
| 33MF | | | | | | | | | |
| 10° | | | | | | | | | |
| 30F | | | | | | | | | |
| 32F | | | | | | | | | |
| 34F | | | | | | | | | |
| 20° | | | | | | | | | |
| 29L | | | | | | | | | |
| 31L | | | | | | | | | |
| 33L | | | | | | | | | |
| 29M | | | | | | | | | |
| 31M | | | | | | | | | |
| 33M | | | | | | | | | |
| ANATOMICAL MEDIUM - 33° | | | | | | | | | |
| 30M | | | | | | | | | |
| 32M | | | | | | | | | |
| 34M | | | | | | | | | |
| ANATOMICAL LONG - 33° | | | | | | | | | |
| 30L | | | | | | | | | |
| 32L | | | | | | | | | |
| 34L | | | | | | | | | |
| Total | | | | | | | | | |

1 x 8 18 LOWER POSTERIOR MOLDS Boxes

| Molds | 200 | 202 | 204 | 206 | 208 | 210 | 212 | 214 | 216 |
|-------------------------|-----|-----|-----|-----|-----|-----|-----|-----|-----|
| 0° | | | | | | | | | |
| 29MF | | | | | | | | | |
| 31MF | | | | | | | | | |
| 33MF | | | | | | | | | |
| 10° | | | | | | | | | |
| 30F | | | | | | | | | |
| 32F | | | | | | | | | |
| 34F | | | | | | | | | |
| 20° | | | | | | | | | |
| 29L | | | | | | | | | |
| 31L | | | | | | | | | |
| 33L | | | | | | | | | |
| 29M | | | | | | | | | |
| 31M | | | | | | | | | |
| 33M | | | | | | | | | |
| ANATOMICAL MEDIUM - 33° | | | | | | | | | |
| 30M | | | | | | | | | |
| 32M | | | | | | | | | |
| 34M | | | | | | | | | |
| ANATOMICAL LONG - 33° | | | | | | | | | |
| 30L | | | | | | | | | |
| 32L | | | | | | | | | |
| 34L | | | | | | | | | |
| Total | | | | | | | | | |

This is a sample of our tooth order form. You may use a photo copy of this form for your first order.

Anterior Mold Comparison

With Dentsply's Bioblend®, Trublend®, Portrait®



084-59-0026/C
 Justi Blend Maxillary
 Centrals Shade Guide
 084-59-0026
 Justi Blend Extended
 Blend Selector in 9 shades



Justi Blend Articulation Chart

| Anteriors | | Posteriors | | | |
|-------------|-----------------------------|----------------|-------------|-----------|-----------------------------|
| Upper Molds | Articulates Lower Anteriors | Zero Degree MF | 10 Degree F | 20 Degree | Anatomical Short, Med. Long |
| A4 | 4 | 31 | 32 | 29L | 30 |
| A7 | 7 | 31 | - | 33M | 32 |
| B3 | 3 | 31 | 32 | 31M | 30 |
| B8 | 8 | 31 | - | 33M | 32 |
| C1 | 1 | 29 | 30 | 29L | 30 |
| C2 | 2 | 29 | 30 | 29L | 30 |
| C5 | 5 | 31 | 32 | 31L | 32 |
| C6 | 6 | 31 | 32 | 31L | 32 |
| D1 | 1 | 29 | 30 | 29M | 30 |
| D2 | 2 | 31 | 32 | 31L | 30 |
| D3 | 3 | 31 | 32 | 29L | 30 |
| D6 | 6 | 33 | 34 | 33L | 34 |
| E1 | 1 | 29 | 30 | 29L | 30 |
| E2 | 2 | 31 | 30 | 29L | 30 |
| E3 | 3 | 31 | 32 | 31L | 30 |
| E6 | 6 | 31 | 32 | 31L | 32 |
| E8 | 8 | 33 | 34 | 33L | 32 |
| F1 | 1 | 29 | 30 | 29L | 30 |
| F4 | 4 | 31 | - | 31M | 30 |
| G5 | 5 | 31 | 32 | 31L | 32 |
| G7 | 7 | 31 | 32 | 31L | 32 |
| H2 | 2 | 31 | 30 | 29L | 30 |
| H4 | 4 | 31 | 32 | 31L | 30 |
| J2 | 2 | 29 | 32 | 29L | 30 |
| J3 | 3 | 29 | 32 | 31L | 30 |
| J5 | 5 | 31 | 34 | 33L | 32 |
| K4 | 4 | 31 | - | 29L | 30 |

| Justi Blend® Bioblend® IPN® Bioblend® Trublend®, Portrait® | | Justi Blend® Bioblend® IPN® Bioblend® Trublend®, Portrait® | |
|---|-----|---|-----|
| Uppers | | Uppers | |
| A4 | 12E | F4 | 43F |
| A7 | 12G | G5 | 45F |
| B3 | 32E | G7 | 45H |
| B8 | 32G | H2 | 55D |
| C1 | 21C | H4 | 55F |
| C2 | 21D | J2 | 62D |
| C5 | 21E | J3 | 55D |
| C6 | 21X | J5 | 62G |
| D1 | 22C | K4 | 75E |
| D2 | 13E | Lowers | |
| D3 | 22E | 1 | C |
| D6 | 22H | 2 | F |
| E1 | 42C | 3 | H |
| E2 | 42D | 4 | N |
| E3 | 42F | 5 | O |
| E6 | 42G | 6 | P |
| E8 | 42H | 7 | R |
| F1 | 43D | 8 | S |

Shade Comparison With Major Brands

| Color | Justi Blend | Bio-blend® | Bio-form® | Tru-Blend™ | Vita® | Portrait™ | Ivoclar Vivadent |
|--------------|-------------|------------|-----------|------------|-------|-----------|------------------|
| White/Gray | 200 | 100 | B59 | T2 | A-1 | P-1 | - |
| Yellow/Pink | 202 | 102 | B53 | T3 | A2-B2 | P-2 | 1A |
| Yellow/Pink | 204 | 104 | B92 | T4 | A2 | P-2 | 1A, 2A |
| Gray/Pink | 206 | 106 | B69 | T5 | C2-D3 | P69 | - |
| Yellow/Brown | 208 | 108 | B62 | T12 | A3-D3 | P-3 | 1C |
| Gray | 210 | 110 | B66/93 | T15 | A35 | P33/P66 | 1C |
| Gray/Brown | 212 | 112 | B67 | T16 | B3 | P67 | - |
| Gray/Brown | 214 | 114 | B83 | T19 | A4-B4 | P3-5 | 3A |
| Gray/Brown | 216 | 116 | B83 | T19 | A-4 | P3-5 | 2B |

This shade comparison chart is to be used in selecting approximate color matches. Slight variations occur due to the difference in the texture and character of the respective materials. This chart should be used as an aid to the skill and preferences of the dentist or technician who must make the final shade choice.

Justi Blend Anteriors match the central, lateral and cuspid, right & left respectively. Each central, lateral and cuspid are the same shade in three different tones thereby creating a natural blended effect.

Justi Blend Posteriors match the right or left cuspid body in tone and shade.

Justi Blend teeth meet ANSI/ADA Specification 15 published in 1999 for Certification.

Instructions for use of synthetic resin teeth

Note: If you wish written instructions for obtaining optimum results with our teeth, contact:

Justi Products, A Brand Division of American Tooth Industries,
 1200 Stellar Drive, Oxnard, CA 93033-2404
 Toll Free (800) 235-4639, California (805) 487-9868
 Fax (805) 483-8482

Distributed by:

MANUFACTURED BY:



Justi Products

A Brand Division of American Tooth Industries

1200 Stellar Drive, Oxnard, CA 93033-2404

(805) 487-9868 • FAX (805) 483-8482

E-mail: info@americantooth.com • Website: www.americantooth.com